



Cleaning your device has never been easier thanks to the water resistant (IP55) feature of Pulmogine®



Silent and portable

Developed to support a better quality of life

Pulmogine®

VIBRATING MESH NEBULIZER


The Missing Piece To Your Inhaled Drug Delivery




The Pulmagine® Mesh Nebulizer sets a novel platform for you drug-device combination product, delivering a wide range of liquid formulations.



No more restrictions for your inhaled medications



Solutions
Salbutamol
2.5ml in 3.6 min*



Suspensions
Budesonide
2ml in 3.5 min*




Biologics
Interferon Beta-1b
1ml in 3.8 min*



High Viscosity Drugs
Liposome
2ml in 5.3 min*

*Data from in-house testing at HCmed Innovations.



Enhance drug delivery performance and achieve optimal results using our customization platform



Customization allows to adjust parameters to provide desired outcome

Smaller MMAD, associated to higher FPF and FPD, ensures greater central and peripheral lung deposition.



Pulmagine® provides affordable top-quality inhalation therapy



Micro Synchronized Delivery Technology™ (MSDT) atomizes a wide range of drugs



Shorter nebulization time and higher delivery rate



Smart design certifies leak-proof feature



Biocompatibility test results guarantee patients' safety

It is easy to use for both adults and children



Pulmagine® Accessories



Adult Mask



Pediatric Mask



Mouthpiece